Q4 CURING LIGHT SPECIALS 2025



BUY AMERICAN

The Best in LED Curing Light Technology



"Excellent Intraoral Access with a Slim 90° head"

- Clinicians Report**









TC-24a LED Corded TC-CL II LED Cordless



TC-01 LED Corded



TC-24 LED Corded

BUY ANY SPRING CURING LIGHT GET THREE 25/PK DIAMONDS



Radiometer (Curing Light Meter) 662

Analog radiometer to check cure lights output from 100-3000 mw/cm2. No battery required.

- Sensitivity Range: 400-500nmIntensity: 0-3000 mW/cm2
- Accuracy: 1 / 50/-



www.springhealthproducts.com

Sold thru your Local Dental Distributor



Must be on one Invoice. Invoices Cannot Be Combined! *No Charge Goods must be of equal or lesser value. Email a copy of you invoice dated 09/01/25-12/31/25 to: customercare@springhealthproducts.com. Include you name, address, email address, phone and choice of free goods shape & grit. Must be received by 01-15-26. MUST INCLUDE CODE "0425".