

# Q2 National Offers 2026

Free goods will ship with order | Offers run 4-1-2026 through 6-30-2026

Use code: **VAQ2\_ENDO**

## NEW! Endocator®

- **Rapid reliable assessment** – Utilizes bioluminescence technology to detect residual contamination invisible to microscopes and radiographs, offering an objective, quantitative score of cleanliness
- **Increased efficiency** – Optimize your disinfection protocol with measurable verification of a thoroughly clean canal
- **Improved confidence** – Provides visible evidence of canal cleanliness, to reduce the likelihood of retreatments

### Endocator® Kit

407330 | 1 x Endocator Device  
1 x USB Wall Charger  
1 x USB Cable  
1 x IFU  
1 x Quick Start Guide

Buy 1 Kit  
Get:

(1) Endocator® Test Swabs (20pk)  
407731  
**\$150 Value (MSRP)**

FOR THIS PRODUCT ONLY - Submit invoice to [orders@vista-dental.com](mailto:orders@vista-dental.com) by 6/17/26 to redeem



## RE-GEN™ Bioactive EndoPutty

- **Ready to use:** A single-component, easy-to-use putty that's highly radiopaque
- **A faster set:** Cures in under an hour, enabling quicker treatments
- **Moldable and versatile:** Firm yet adaptable to fit various clinical needs and anatomies

RE-GEN™ Bioactive EndoPutty

504756 | 1 x 1g prefilled syringe

Buy 1  
Get 1  
**FREE**



## RE-GEN™ Bioactive EndoSealer

- Ready to use: Injectable, single-component, bioactive paste with smooth flow handling
- A faster set: Fast and efficient set time of under 1 hour
- Dimensionally stable: Low solubility with high wash out resistance

RE-GEN™ Bioactive EndoSealer

504753 | • 1 x 2g prefilled syringe  
• 1 x 20 Cup MST Visco-Tips™

Buy 1, Get 1  
MST Visco-Tip™  
Polymer Tips  
(50 cup) 312300  
**FREE**



# Q2 National Offers 2026

Free goods will ship with order | Offers run 4-1-2026 through 6-30-2026

Use code: **VAQ2\_ENDO**

## Endo Sealer Solvent™

The first endodontic irrigant designed for bioceramic sealer retreatment and calcium hydroxide removal.

- **Targeted Efficacy:** Engineered to selectively remove sealer and calcium hydroxide while being gentle on dentin
- **Versatile:** Suitable for all bioceramic sealers and calcium hydroxide
- **Time-Efficient:** Streamlines two visit endodontic procedures & retreatment procedures

### Endo Sealer Solvent™ Bottle

317033 | 1 x 8oz (240mL) Bottle



Buy 1 Kit  
Get:

(1) Appli-Vac™ 1" 30GA  
Irrigating Tips (100pk)

315530



## EDTA

### 17% Solution

Engineered with lower pH for optimal performance—pH matters more than concentration

- Opens dentin tubules
- Removes smear layer and dentin mud
- Prepares dentin walls for better ADHESION of filling materials

### EDTA 17% Solution Bottle

317011 | 1 x 16oz (480mL) Bottle



Buy 1  
Get

6% Sodium  
Hypochlorite  
16oz Bottle  
317008



## SmearOFF™

### 2-in-1 EDTA/CHX Mix

- EDTA / CHX mix formulated with wetting agents and surface modifiers
- Will NOT form a precipitate when used with NaOCl
- Does not require a rinse

### SmearOFF™ 16oz Bottle

317111 | • 1 x 16oz Bottle



Buy 1  
Get

Chlor-XTRA™  
16oz Bottle  
503800



## Appli-Vac™ 1" Irrigating Tip

- Luer Lock Hub
- Slotted and side-vented for safe irrigation

### Appli-Vac™ 1" Irrigating Tip

- 315523 Appli-Vac™ 1" Irrigating Tip 23ga 100pk 100 x 23ga (Blue)
- 315527 Appli-Vac™ 1" Irrigating Tip 27ga 100pk 100 x 27ga (Yellow)
- 315530 Appli-Vac™ 1" Irrigating Tip 30ga 100pk 100 x 30ga (Purple)

Buy 3  
Get 1  
FREE

